

## Our Experience – Bed Slide Sheets



Many companies manufacture Bed Slide Sheets, and they are one of the original Manual Handling products of the last 30 years. They come in a variety of colours and materials and are generally made from a nylon, parachute type material. Different manufacturer's materials have different characteristics.

The Pelican range has four colours and two sizes:

Colour	Code	Comments
Red	507Q, 507Q2	Quiet, breathable material that can be left under patients.
Pink	507S, 507S2	Plastic bag type noise. Not breathable. The best for slipperiness.
Pale Blue	507QB, 507QB2	Quiet, breathable material that can be left under patients.
Royal (Dark) Blue	507RB, 507RB2	Slightly noisy material. Not breathable. Should not be left under patients.

We also have two sizes: **1m x 1.5m** or **2m x 1.5m**. Bed Slide Sheets are normally **folded in half**. Often the smaller size can end up being too small for some people, so we recommend purchasing the larger sheets. The sheets are inserted underneath someone in bed by the log-rolling Manual Handling technique. With regard to weight limits, this is an area where the ability of the Carer is what matters. The strength of the material is almost irrelevant because:



- 1. The sheet will never be taking someone's full weight as it is only used to help turn the patient**
- 2. The material's strength will far outweigh the strength of the Carer**

A few items to be aware of are:

### **Bed Slide Sheet Holder** Item # 509H, 509HL

The original smaller size will hold two 1.5m x 1m sheets, or one 1.5m x 2m sheet. The larger size is designed to hold two 1.5m x 2m sheets. They are made from a soft clear vinyl material and can be wiped clean or hand or machine washed. They can be hung from a hook, door knob, cupboard rail or the bed frame. It is good to keep the Bed Slide Sheets in a dedicated place to help prevent people slipping on them if they fall on the floor, and also to help prevent cross-infection.



### **Bed Slide Sheet Gripper** Item # 507BSSG

This product is a long piece of firm foam which you can wrap the outer edge of the Bed Slide Sheet around. This makes it easier for the Carer to hold onto, relieving the pressure on the Carer's fingers.



**Fingernails:** In the past we have been made aware that some Carer's fingernails can 'spike' the material. This is something to be aware of.

**Static Electricity:** The question of static electricity has also been raised with using Bed Slide Sheets. Testing showed that the 507S pink Bed Slide Sheet material created 1,200 Volts. This was done by rubbing it together and measured using a fully calibrated ACL Digital Static Field Analyser. Things like plastic trays, hoses, rubber fittings, barrier gowns and drapes, polycotton gowns and drapes can generate up to 20,000 – 30,000 Volts. Normally to create a spark two pieces of metal need to rub together. Staff who have had static shocks have been wearing polyester jumpers and working on a type of carpet that has contributed to the static.

**Cleaning:** Wash and do not use fabric softener. Re-spray with Helmar Silicone Spray to increase the slipperiness.

### **Slide & Turn Sheet, Item # 509**

The main advantage is that this can be left under a patient during the night, because it is made from a breathable material which will not increase sweating, and the Carers do not need to log-roll the patient each time they are turned, which is what they need to do with normal Bed Slide Sheets. It has a black breathable foam pad which anchors it to the bed to help it stay in position.



Also available: **Slide & Turn Hoist Sheet** designed to allow a hoist to do the heavy lifting for heavy patients.

### **Key points:**

- **4 colours** with different features
- Bed Slide Sheets **do not take someone's full weight**
- Weight limit is **limited to the Carer's ability and strength** and not to the strength of the material
- Breathable **Slide & Turn Sheet** can be left under patients
- Breathable **Slide & Turn Hoist Sheet** is available for heavier patients